### P-23 Low Energy Neutron Team

David Bowman
Vincent Yuan
Seppo Penttila
Scott Wilburn
Greg Mitchell
Gil Peralta
Shelley Page
(visitor)

#### **About the NPDGamma Experiment**

$$\overrightarrow{\mathbf{n}} + \mathbf{p} \to \mathbf{d} + \gamma$$

#### The NPDGamma Experiment at





**LANSCE** is the only facility worldwide that can host this nuclear physics experiment.

NPDGamma is under construction and will begin data collection in 2003.

# Measurement of the Parity-Violating Gamma Asymmetry $A_{\gamma}$ in the Capture of Polarized Cold Neutrons by Para-Hydrogen, $\vec{n}+p \rightarrow d+\gamma$

J.D. Bowman (Spokesperson), G.L. Greene, G.S. Mitchell, S.I. Penttilä, W.S. Wilburn, V.W. Yuan Los Alamos National Laboratory

T.A. Case, S.J. Freedman *University of California, Berkeley* 

T.B. Smith *University of Dayton* 

E.I. Sharapov

Joint Institute for Nuclear Research, Dubna

G.L. Jones *Hamilton College* 

M. Gericke, I. Kuznetsov, H. Nann, W.M. Snow *Indiana University* 

S. Ishimoto, Y. Masuda, K. Morimoto KEK National Laboratory, Japan

S.A. Page, W.D. Ramsay University of Manitoba and TRIUMF

T.E. Chupp, K.P. Coulter *University of Michigan* 

T.R. Gentile

National Institute of Standards and Technology

F.W. Hersman, M.B. Leuschner *University of New Hampshire* 

R.D. Carlini

Thomas Jefferson National Accelerator Facility

http://p23.lanl.gov/len/npdg/

#### What is NPDGamma measuring?

The laws of physics describe four interactions:
 gravity
 electromagnetism
 strong force
 weak force

The strong force binds hadrons together, and is the primary interaction between protons and neutrons.



The weak interaction between hadrons

(i.e. protons and neutrons)

is not well-measured or understood.

How do we study the hadronic weak interaction?


To see its sometimes very small effects, a powerful technique is to use parity violation.

force	strong	weak	
strength	$\sim 10^6$ 1		
parity?	conserved	violated	

(Necessary to observe  $\Delta s = 0$ , nonleptonic weak interactions.)

$$\overrightarrow{n} + p \rightarrow d + \gamma$$
 (2.2 MeV)

NPDGamma will measure  $A_{\gamma}$ , the parity-violating asymmetry in the distribution of emitted  $\gamma$ 's



mirror image

If the up/down  $\gamma$  rates differ, parity is violated (PV  $\longrightarrow$  signature of the weak interaction)

Expected asymmetry  $\approx 5 \times 10^{-8}$ 

Goal experimental error:  $0.5 \times 10^{-8}$ 

Range of  $Z, W^+, W^-$  bosons is 0.002 fm


But nucleon interactions take place on a scale of 1 fm (short range repulsion)




the weak force interaction between nucleons and hadrons is a meson exchange

At low energies (< 300 MeV) mesons are the appropriate degree of freedom

Meson exchange model is a successful picture of strong interactions between nucleons (describes to a few % n-p/p-p scattering cross-sections)



 $A_{\gamma}$  is a clean measurement of  $H_{\pi}^{1}$ :  $A_{\gamma} \approx -0.045 H_{\pi}^{1}$ 



NPDGamma will provide a measurement with improved statistical precision compared to  $^{18}$ F results, with no uncertainties from many-body calculations or nuclear structure effects.


#### What does NPDGamma need?







#### **NPDGamma Experimental Setup**




NPDGamma is a funded experiment (\$4.8M)

#### DOE & NSF

also Hamilton College, Indiana University, KEK, University of Manitoba, University of Michigan, NIST, University of New Hampshire, TJNAF

#### **Overview**



Surface Profile Analysis Reflectometer (SPEAR) is used with an unpolarized neutron beam to study solid/solid, solid/liquid, solid/gas, and liquid/gas interfaces.

Greg Smith, 505-665-2842, gsmith@lanl.gov and Jaroslaw Majewski, 505-667-8840, jarek@lanl.gov

Moderators

1, 2, 16: high-resolution H<sub>2</sub>O decoupled

- Low-Q Diffractometer (LQD) is designed to study structures with dimensions in the range from 10 to 1000 Å. It measures a broad Q-range in a single experiment without physical changes to the instrument. Rex Hjelm, 505-665-2372, hjelm@lanl.gov
- AS ERIX will provide a polarized neutron beam for studies of magnetic materials, using reflectometry and diffraction, and includes application of high magnetic fields

Mike Fitzsimmons, 505-665-4045, fitz@lanl.gov

- FP12 will be used for a fundamental nuclear physics experiment to precisely measure the asymmetry of the emission of gamma rays from the capture of polarized neutrons by protons. David Bowman, 505-667-7633, bowman@lanl.gov
- IN500 is a protoype instrument under development employing novel techniques to enhance inelastic coldneutron spectroscopy at spallation neutron sources. Margarita Russina, 505-667-8841, russina@lanl.gov Ferenc Mezei, 505-667-7633, mezei@lanl.gov
- **Detector for Advanced Neutron Capture Experiments** (DANCE) will be used for the study of neutron capture on radioactive nuclei in support of the stockpile stewardship program and for nuclear astrophysics. John Ullmann, 505-667-2517, ullmann@lanl.gov
- Protein Crystallography Station (PCS) is a single-crystal FP15 diffractometer designed for structure determinations of large biological molecules. Paul Langan, 505-665-8125, langan\_paul@lanl.gov Benno Schoenborn, 505-665-2033, schoenborn@lanl.gov
  - **PHAROS** is a high-resolution chopper spectrometer designed for studies of Brillouin scattering, magnetic excitations, phonon densities of state, crystal-field levels, and chemical spectroscopy and measurements of  $S(Q,\omega)$ . Robert McQueeney, 505-665-0841, mcqueeney@lanl.gov

Neutron Powder Diffractometer (NPD) allows for studies of complex structures, internal strain measurements, and phase transformation. Don Brown, 505-667-7904, dbrown@lanl.gov

FP2 Spectrometer for Materials Research at emperature and Stress (SMAR S) will allow measurements of spatially resolved strain-fields, phase deformation and load transfer in composites, the evolution of stress during temperature (or pressure) fabrication, and the development of strain during reactions (such as reduction, oxidation, or other phase transformations).


Mark Bourke, 505-665-1386, bourke@lanl.gov

- High Intensity Powder Diffractometer (HIPD) is designed to study the atomic structure of materials that are available only in polycrystalline or noncrystalline forms.
  - Robert Von Dreele, 505-667-3630, vondreele@lanl.gov
- High-Pressure-Preferred Orientation (HIPPO) instrument is a new high-intensity powder diffractometer for high-pressure and texture

Kristin Bennett, 505-665-4047, bennett@lanl.gov and Robert VonDreele, 505-667-3630, vondreele@lanl.gov

- FP5 is used to study the Doppler shift and broadening of low-energy nuclear resonances in materials under extreme conditions and for structural studies using transmission Bragg diffraction. Vincent Yuan, 505-667-3939, vyuan@lanl.gov
- Single Crystal Diffractometer (SCD) has been used to study the structure of organometallic molecules, unique binding of H2 crystal structure changes at solid-solid-phase transitions, magnetic spin structures, twinned or multiple crystals, and texture. Yusheng Zhao, 505-667-3886, yzhao@lanl.gov
- Filter Difference Spectrometer (FDS) is designed to determine energy transferred to vibrational modes in a sample by measuring the changes in the energies of the scattered neutrons. Juergen Eckert, 505-665-2374, juergen@lanl.gov

# NPDGamma building FP12 to be ready for: commissioning run Fall 2002 production data taking 2003



(MOU between LANSCE and P Divisions)

### Flight Path 12 Construction Progress

in-pile guide



shutter



frame overlap chopper




### <sup>3</sup>He Spin Filter

Optical pumping of Rb vapor, which polarizes <sup>3</sup>He by hyperfine spin-exchange collisions.


Neutron beam is polarized by passing through the cell. Antiparallel spin neutrons are absorbed.




Fall 2000 Test Run:  $^3$ He polarization of 26.5%  $\longrightarrow n$  polarization of 30-70% for 2-10 meV



## Liquid para-hydrogen target 20 $\ell$ , Mg-Al cryostat window, <sup>6</sup>Li liner




n cross-sections: ortho-  $(\uparrow\uparrow)$  and para-  $(\downarrow\uparrow)$  hydrogen  $\Box$  ortho- scattering,  $\circ$  para- scattering, -np capture





at 17K, ortho- fraction is 0.03%

Para-hydrogen: necessary to allow neutron capture and to preserve neutron polarization upon scattering


## CsI(Tl) and Photodiode $\gamma$ Detectors 48 of these detectors will be used in the full experiment







## PV Asymmetry Measurements on Cl, La, Cd Test Run Fall 2000



#### Raw PV Asymmetries $(\times 10^{-6})$

<sup>35</sup> CI	$-7.68 \pm 2.17$
<sup>139</sup> La	$-5.88 \pm 2.35$
<sup>113</sup> Cd	$+1.94 \pm 1.48$

#### Physics Asymmetries $A_{\gamma}$ (×10<sup>-6</sup>)

	<sup>35</sup> Cl	<sup>113</sup> Cd	<sup>139</sup> La
Leningrad	$-27.8 \pm 4.9$	$-1.3 \pm 1.4$	$-17.8 \pm 2.2$
ILL	$-21.2\pm1.7$	-	-
LANSCE	$-23.1 \pm 6.5$	$+5.8 \pm 4.4$	$-17.1 \pm 6.8$

(LANSCE results to be published.)

#### **NPDGamma Status**

- FP12 flight path is under construction.
- Experiment is under construction.
   10% scale apparatus tested Fall 2000.
   Alignment scheme & monitors tested Fall 2001.
   All crucial components demonstrated.
- Test runs indicate design is sufficient for target  $A_{\gamma}$  experimental error,  $0.5 \times 10^{-8}$ .
- Potential systematic errors studied extensively.
- NPDGamma will make a clean measurement of  $H_{\pi}^{1}$ , the most fundamental weak N-N coupling.

#### **NPDGamma Schedule**

January 2002 Start beamline installation

November 2002 Commission beamline through ER1

January 2003 Installation of guide complete

June 2003 FP12 Commissioning Run

August 2003 Commission entire experiment

Fall 2003 Begin data taking

December 2003 Data match existing stat. precision on  $H^1_\pi$ 

2004 / 2005 Continue data taking for 10% measurement